# MKG SPEZIAL

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
All-over length of instruments may vary depending on length of working part and neck!

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„Saw-toothing“ with 4 blades for fast cutting of bone substance, **ZrN-coating as wear-protection coating**

- high corrosion resistance
- easy reprocessing/ cleaning
- very smooth surface ensures a constant cutting capacity during the intervention
- reduced wear ensures more efficiency

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Particularly smooth-running, fast cutting of bone substance and tooth structure, **ZrN-coating as wear-protection coating**

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lengthened working part for a wider cutting surface

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Surgical spiral toothing with cross cut

|                             |                             |                             |                             |                             |
|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|
| 506 <b>104</b> 408 297 016T | 506 <b>205</b> 408 297 016T | 506 <b>316</b> 408 297 016T | 506 <b>104</b> 409 297 021T | 506 <b>205</b> 409 297 021T |
| 1,6 / 9,2                   |                             |                             | 2,1 / 9,2                   |                             |

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Fine 10-blade ball cutter, especially smooth-running

|                            |                            |                     |                            |
|----------------------------|----------------------------|---------------------|----------------------------|
| 500 <b>104</b> 001 251 027 | 500 <b>205</b> 001 251 027 | 500 104 001 251 031 | 500 <b>205</b> 001 251 031 |
| 2,7 / 2,7                  |                            | 3,1 / 3,1           |                            |

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Fine surgical cross-cut, particularly smooth-running

|                     |                     |                            |                     |                     |                     |                     |
|---------------------|---------------------|----------------------------|---------------------|---------------------|---------------------|---------------------|
| 500 104 001 258 018 | 500 104 001 258 021 | 500 <b>104</b> 001 258 023 | 500 104 001 258 027 | 500 104 001 258 031 | 500 104 001 258 035 | 500 104 001 258 040 |
| 1,8 / 1,8           | 2,1 / 2,1           | 2,3 / 2,3                  | 2,7 / 2,7           | 3,1 / 3,1           | 3,5 / 3,5           | 4,0 / 4,0           |

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HPL TC-Cutter · surgery diamonds

**258**

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Surgical cross-cut, particularly smooth-running

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Pilot drill are used for the initial determination of the drilling depth. Afterwards the reaming of the drilling of the implant site is made by the final drill, for an optimal fit of the implant body.

- stainless special steel with a low cracking tendency
- sharp 2-blades-construction with big flutes for an easy chip removal
- s-shaped blade tip for easy penetrating
- relief-ground flutes for a low friction on the bone
- **laser-depth marking for accurate drilling (after each 2 mm), 2-blades after 8 mm, 3-blades after 6 mm drilling depth**

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

**Optimum speed:** 800 min<sup>-1</sup> – 1.200 min<sup>-1</sup>. Do not exceed the maximum of 2N.  
Always use with spray, 50ml/min.

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Due to insufficient oral hygiene might result a plaque formation on the implants and cause inflammations. Finally, this osteolysis causes a damage of the implant bed. The implant body is partially exposed.

The accurate peri-implantitis finishers are outstandingly suitable to smooth the implant surfaces and support the effort to preserve the implant.

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For removal of the bone cover during the root tip resection.  
By the maximum drilling depth of 5.0 mm and the reduced total length, those instruments are also suitable for the posterior region (premolars/molars)

|      | 330 204 227 040 | 330 204 227 050 | 330 204 227 060 | 330 204 227 070 |
|------|-----------------|-----------------|-----------------|-----------------|
| L    | 5,0             | 5,0             | 5,0             | 5,0             |
| Ø mm | 4,0 / 2,9       | 5,0 / 3,9       | 6,0 / 4,9       | 7,0 / 5,9       |

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For the preparation of bone cylinders for bone harvesting.  
By the maximum drilling depth of 8 mm and the reduced total length, these instruments allow a better visibility on the worksurface and reduce consequently the risk of injury

|      | 330 204 228 040 | 330 204 228 050 | 330 204 228 060 | 330 204 228 070 |
|------|-----------------|-----------------|-----------------|-----------------|
| L    | 8,0             | 8,0             | 8,0             | 8,0             |
| Ø mm | 4,0 / 2,9       | 5,0 / 3,9       | 6,0 / 4,9       | 7,0 / 5,9       |

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For explanting.

- the maximum drilling depth allows a safe explanting of the commercially available implants up to the size of the diameter of 4.8 mm and a maximum total length of 18.0 mm
- the milled down windows allow a good visibility on the implant body and support the cooling for an atraumatic explantation

|      | 330 204 229 035 | 330 204 229 050 | 330 204 229 060 |
|------|-----------------|-----------------|-----------------|
| L    | 18,0            | 18,0            | 18,0            |
| Ø mm | 3,5 / 2,9       | 5,0 / 3,9       | 6,0 / 4,9       |

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**Recommendations for Trepan bur:** To use with contra-angle handpiece (minimum speed reduction 10:1) and permanent, sterile cooling. Optimum Speed: 800–1.000 min<sup>-1</sup> Max. Speed: 5.000–6.000 min<sup>-1</sup>

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Crown cutter RR for all metallic materials

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Crown cutter Larry Black for all types of ceramic crowns

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Separates all ceramic materials, easily and quickly

- conical, extended working part
- interrupted cut provides a quick removal
- improved cooling capability due to a corrugated profile
- more diamnd grains by an enlarged surface
- high performance diamond coating for max. durability

|                     |  |
|---------------------|--|
|                     |  |
| 800 314 473 528 016 |  |
| 1,6 / 8,0           |  |

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acurata instruments are intended for the dental medical field and must only be used by dentists or respective experts who, due their training and experience, are familiar with the use of the products. The instruments are supplied non-sterile and have to be prepared prior to first use.

## INFORMATION FOR USE FOR DENTAL ROTARY INSTRUMENTS

Use only technical and hygienic correct drives with a stable ball bearing! Carefully clamp the instruments as deep as possible and assure that the desired speed is reached prior to usage. Please note the maximum allowed engine speed as stated on the packaging label. Avoid levering or canting, otherwise there is risk of fracture. Pay attention that there is sufficient cooling (50 ml/min) otherwise there is risk of thermal damage. FG-instruments with a total length > 19 mm and instruments with a diameter of the work element > 2 mm (ISO 020) need extra cooling. Separate blunt, twisted or damaged instruments immediately and do not use them again. Improper usage increases risk, higher wear and inferior work results!

**Contact pressures:** The contact pressures should be light, do not exceed the maximum of 2N. Excessive contact pressures can cause damage to the instruments. In extreme cases the instrument can fracture! At the same time a higher heat development is generated. Overheating may harm the dental pulp, instruments wear faster and surfaces get coarser.

**Remarks:** Polishers and instruments with long working part and respectively long, slender necks (e.g. instruments for micro preparation) or for specific application areas (e.g. crown separator, long surgical cutter): exceeding the maximum allowed speeds can cause oscillations (vibrating) or displacement (whippings). This can result in fracture or heavy damage. Non observance causes an increased security risk! Therefore pay attention to the maximum speed printed on the packaging labels. The respective maximum speed is a theoretical value given on the basis of the geometry of the instruments. The choice of the working speed is determined by the material subsequently to be used, the specific indication, the type of handpiece used and the contact pressures and is at user's discretion. Please ensure a dry, dust-free and sun-protected storage at appropriate room temperature.

### Maximum speeds – table of reference values

**FG instruments (ISO 314)** ISO Ø: 005–016 max. 450.000 min<sup>-1</sup> • 018–021 max. 300.000 min<sup>-1</sup> • 023–031 max. 160.000 min<sup>-1</sup> • 033–055 max. 120.000 min<sup>-1</sup>

**HP and RA instruments (ISO 104 / 204)** ISO Ø: 003–023 max. 50.000 min<sup>-1</sup> • 025–040 max. 40.000 min<sup>-1</sup> • 045–060 max. 30.000 min<sup>-1</sup> • 060–070 max. 25.000 min<sup>-1</sup> • 070–250 max. 25.000 min<sup>-1</sup> • 300–380 max. 20.000 min<sup>-1</sup> • 450 max. 15.000 min<sup>-1</sup>

**Polishing instruments** all ISO Ø: max. 5.000–10.000 min<sup>-1</sup>

The information on the label on the packaging applies. Optimum speeds for all instruments depend on the material to be processed. As a general rule 40 – 50% of the maximum instruction is recommended!

## RECOMMENDATIONS FOR USE FOR ROOT CANAL INSTRUMENTS

K-Reamer: push and rotary motion clockwise up to max. 90°. K-File: filing motion clockwise up to max. 45°, Hedstroem: push and pull (with or without rotary) motion up to max. 90°. Usage in ascending order of size (do not skip any size), High flexible Nickel-Titanium instruments are especially for heavy curved root canals. The file has not to be bent before use, Protect the instruments against aspiration (rubber dam – safety chain – safety cord), Silicone stoppers have to be removed prior to reprocessing root canal instruments, Clean the instruments properly before sterilization, Sterilization in autoclave may corrode steel instruments, Do not use bent or damaged instruments.

**Cleaning, disinfection, sterilization and storage of medical products – refer to our Recommended Hygiene Procedures, also available on our website:**  
[www.acurata.de/en-US/products/hygiene-recommendations](http://www.acurata.de/en-US/products/hygiene-recommendations).

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# RECOMMENDED HYGIENE PROCEDURES FOR DENTISTRY

rev. 07/15

Manufacturer's information on the reprocessing of medical products semi critical B and critical B according to the RKI recommendation "Prevention of infections in dentistry hygiene requirements", the KRINKO guideline "Requirements on the hygiene for reprocessing of medical products" and the DIN EN ISO 17664.

**Manufacturer:** acurata GmbH & Co. KG · Schulstraße 25 · 94169 Thurmansbang · Tel.: +49(0)8504-91 17 0 · Fax.: +49(0)8504-91 17 90 · info@acurata.de · www.acurata.de

**Products:** The present manufacturer's information applies to all instruments supplied by acurata GmbH & Co. KG that are used for surgical, periodontal or endodontic treatments. These are rotating tungsten carbide, diamond and polishing instruments as well as oscillating instruments made of stainless steel or nickel-titanium. The product range includes exclusively instruments delivered in non-sterile condition. These have to be prepared prior to first use (begin with step 1) and prior to any further use (begin with step 2).

**Limited number of reprocessing cycles:** The end of a product's service life is determined by damage and wear caused by use. Frequent reprocessing does not affect the performance of these instruments. Disposable products (marked  on the label) shall not be reused or reprocessed.

**Basic note:** Observe the legal provisions concerning the reprocessing of medical products valid in your country (e.g. in Germany www.rki.de). The manufacturer assures that the reprocessing methods detailed below are appropriate for the reprocessing of the mentioned groups of instruments for their reutilization in regard of their intended use. The operator is responsible that the applied methods of reprocessing, based on his risk assessment with the used equipment, material, process parameters and staff achieve the required results for the intended use. To guarantee this, routine controls of the validated mechanical and/or manual reprocessing methods are necessary. Any deviation from the validated method below detailed must be checked and released by the operator to ensure effectiveness and to avoid possible adverse consequences. **The hygiene procedure recommendation is available on our website in its current version:**  
<https://www.acurata.de/en-US/products/hygiene-recommendations>

## 1 Preparation incl. storage and transportation

For first use preparation begin with step 2. Place instruments immediately or at the latest one hour after use on a patient, in a cleaning/disinfection tank filled with a suitable detergent/disinfectant (non-fixing/ aldehyde-free, e.g. BIB forte eco). Set up the drill bath acc. to manufacturer's instructions; for BIB forte eco mix concentrate with water - first the water, then add concentrate. Cover the tank. Pay attention to the application time (e.g. BIB forte eco 0.5% 60 min.). The transport of the instruments to the place of preparation should be made in a contamination protected cleaning/disinfection tank.

## 2 Cleaning and disinfection

According to the directive of the Robert Koch Institute (RKI) and Commission for Hospital Hygiene and Infection Prevention (KRINKO) it is preferential that the preparation of semi-critical B products is carried out mechanically; critical B products shall in any case be prepared mechanically. For products with long, tight lumina or cavities the cleaning has to be performed mechanical. For root canal instruments silicone stoppers have to be removed prior to reprocessing.

### Mechanical cleaning – validated method

**Equipment:** Cleaning brush, mechanical washer/disinfector acc. to EN ISO 15883 (RDG) (e.g. Miele with Vario TD program), detergent (e.g. 0.5% cleaner Neodisher mediclean), acurata bur stand made of stainless steel.

**Method:** Remove instruments from instrument stand / cleaning-disinfection-tank immediately before mechanical reprocessing and brush off all visible contamination in a cold water bath. Place the instruments into the opened bur stand. Start mechanical cleaning according to the instructions of the manufacturer of the machinery and of the detergent. The following process is validated: Program Vario TD: 2 min. precleaning, 5 min cleaning at 55 °C with detergent, 3 min. neutralizing, 2 min. intermediate rinsing, final rinsing with appropriate (VE-) water 5 min. at >90 °C.

### Thermal disinfection in a validated washer/disinfector

Perform the mechanical cleaning (see above, e.g. Miele washer/disinfector with Vario TD program) incl. thermal disinfection with the instruments fixed in a bur holder. The manufacturer's instructions for the device must be observed. For validated washer/ disinfectors the disinfection is demonstrably assured. acurata products are thermostable up to 134 °C.

### Manual cleaning and disinfection – standardized method

**Equipment:** Cleaning brush (e.g. synthetic brush, sterilizable), ultrasonic bath, detergent & disinfectant with approved efficiency for dental instruments (e.g. BIB forte eco, Alpro Medical), bur stand for rotating oscillating instruments (e.g. acurata bur stand made of stainless steel); the manufacturers' instructions must be observed.

**Method:** Remove instruments from bur stand / cleaning-disinfection tank immediately before manual cleaning. Brush off all visible contamination in a cold water bath. Rinse the instrument and the bur stand under running water. Put the instrument into a suitable strainer element and place it into the ultrasonic unit filled with detergent & disinfectant. Perform cleaning and disinfection according to the instructions of the manufacturers of the ultrasonic bath and the detergent and disinfectant; e.g. BIB forte eco 3% - 10 min. at 55 °C tested to EN 14476. After the application time rinse the instrument thoroughly with appropriate water (e.g. VE-water). Dry the instruments preferably with medical compressed-air. According to KRINKO the manual reprocessing is finalized by a thermal disinfection in a steam sterilizer. Follow manufacturer specifications.

**Visual examination** with a suitable enlarger to ensure that the instrument is clean and undamaged (an enlargement of 8x -10x is recommended). If after reprocessing still residues of contamination are visible, repeat the cleaning and disinfecting process until no visible contamination is left. Instruments showing defects are to be discarded immediately (e.g. missing diamond coating, blunt and chipped blades, deformations, corroded surfaces or non-removable residual contamination).

## 3 Final reprocessing steps – sterile packaging and sterilization – Medical products critical B – validated method with moist heat:

**Equipment:** Steam sterilizer Co. MMM Selectomat HP acurata bur stand made of stainless steel, transparent sterilization bag (Steriking o. VP Stericlin), sealed seam device Co. Hawo

**Packaging:** Prior to sterilization place the instruments in the bur holder and pack them altogether doubly in sterilization bags and weld them with a sealed seam device. The instruments must be protected. For the packing an appropriate standardized method has to be applied.

**Sterilization:** An effective steam sterilization of the packaged instruments is proven successfully in the pre-vacuum steam sterilization method with the following minimal parameters: 3 pre-vacuum phases, 132 °C sterilization temperature, holding time 3 min. (full cycle), drying time 10min. Follow the instructions of the device manufacturer. Note: The products are not suitable for a sterilization in a hot-air sterilizer or chemiclave.

## 4 Transport and storage During reprocessing, the reprocessed products must be protected from recontamination. Further the packed sterile goods must be protected also from dust and moisture.

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## HIGH PERFORMANCE DENTAL INSTRUMENTS

to the highest technological level.

Dentists and dental laboratories now value our practically-oriented, needs-oriented products all over the world. We also listen and ask questions - at the same level as our customers. It is in this context that the highly focused, dedicated perfectionists at acurata manufacture absolutely reliable precision instruments which enable ultra-accurate work in the laboratory and in the dental practice. Reliability in the context of a sensitive collaboration between all of our talents – at the technical and human level.

WWW.ACURATA.DE

📍 acurata GmbH & Co. KG · Schulstraße 25 · 94169 Thurmansbang 📞 Telefon +49-(0) 85 04 - 91 17 0 📠 Fax +49-(0) 85 04 - 91 17 90

# acurata