PRAXIS- INSTRUMENTE

Rotating instruments for Dentistry

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Precision Precisely sized cutting geometries save tooth substance and materials.

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Smooth running A proved smooth running of our instruments relieves the patient during the treatment and save the tooth substance.

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Durability Professional qualification as well as technical ability ensure an extra durability of our instruments.

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
stands for high performance carbide instruments from acurata, Germany. Manufactured from micro-grain tungsten carbide in state-of-the-art CNC machines, these instruments meet the highest demands of dentists and technicians. Precisely milled blade geometries, high efficiency and vibration-free concentricity ensure smooth cutting of dental hard tissues and a long service life. Quality is managed through all manufacturing stages according to DIN EN ISO 13485. This results in atraumatic preparations. The excellent true running of the instruments has a soothing effect on the patient.

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Round bur for the excavating of carious dentin.

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is the further development of the successful Blue Cut TC-bur

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Particularly long neck,
optimum view of the working surface

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
Carbide round burs with zircon nitrite wear protection layer

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- with higher hardness and corrosion resistance
- minimises wear (demonstrably from medical field, e.g. with prostheses, implants or artificial joints)
- diamond-cut precision blades
- nick-free and gentle to enamel structure
- spiral angles ensure low vibration
- spade-shaped face blade with wide flute spaces
- high penetration capacity
- gentle to the pulp due to low working pressure required

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All instruments are made from one solid tungsten carbide work-piece. There are no soldering joints. Under load, the flow of stress is passed on from the instrument head to the shank and subsequently absorbed by the collet chuck. In this way fracture risk is reduced significantly.

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- innovative cutting edge geometries
- cuts even hardest alloys quickly and efficiently
- high concentricity
- one-piece-tungsten carbide construction
- **tip: excellent for crowns of non-precious-alloys**

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Separates all ceramic materials, easily and quickly

- conical, extended working part
- interrupted cut provides a quick removal
- improved cooling capability due to a corrugated profile
- more diamond grains by an enlarged surface
- high performance diamond coating for max. durability

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Crown cutter, FG-grinder for the treatment of ceramic crowns

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For all-ceramic crowns made of ZrO

- slender, extended working part for extensive finishing
- smooth, controlled running
- long lasting high performance diamond coating

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For all-ceramic crowns made of ZrO

- conical, very strong working part for separation cuts
- smooth, controlled running
- long lasting high performance diamond coating

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For ceramic veneering and low melting ceramics

- slender working part
- high cutting performance “super coarse”
- very efficient

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Quickly shaped!

HM-Former for effective treatment of metals, like gold alloys, non-precious alloys or titan directly on the patient. The easy cutting cross-cut enables an efficient moulding or reduction of implants, individual or assembled pivot teeth. With only one instrument you achieve a high removal performance and plane, easily polishable surfaces. The short, grainy chippings are easy to extract. Forming and finishing of the crown margin without changing the instrument for accurately fitting preparation edges and a dense margin fit.

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Optimum speed: 80.000 min⁻¹ – 100.000 min⁻¹.
 Do not exceed the maximum of 2N. Always use with spray, 50ml/min.

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For perfect occlusion. For natural function. For composites, compomers, glass ionomers, metals (including titanium)

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Due to insufficient oral hygiene might result a plaque formation on the implants and cause inflammations. Finally, this osteolysis causes a damage of the implant bed. The implant body is partially exposed. The acurata peri-implantitis finishers are outstandingly suitable to smooth the implant surfaces and support the effort to preserve the implant.

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- with fine cross-cut for composite
- readily identifiable by the orange ring
- optimum speed 20.000–80.000 min⁻¹
with spray (50 ml/min.)

We recommend to use for the following polishing our diamond interspersed composite polishers Composoft CS (p. 62).

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After an orthodontic treatment, the brackets have to be removed. The frequently left adhesive residues have to be removed of the enamel. Therefore the HOS 016 have been specially developed. They remove the adhesive residues smoothly and rapidly and save the enamel. The stepless end of the bur excludes scratches and grooves on the enamel. The best tungsten carbide and a precise manufacturing ensure a smooth running and profitability.

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Characteristics:

- fast cutting version
- smooth running
- very efficient
- one-piece-tungsten carbide construction

Characteristics:

- fast cutting version
- smooth running
- very efficient
- one-piece-tungsten carbide construction

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Adhesive Remover, easy removal of adhesive residues –
ZrN-coating as wear-protection coating

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acurata surgery tungsten carbide instruments are primarily characterized by fast separation and they are particularly smooth running.

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„Saw-toothing“ with 4 blades for fast cutting of bone substance,

ZrN-coating as wear-protection coating

- high corrosion resistance
- easy reprocessing/ cleaning
- very smooth surface ensures a constant cutting capacity during the intervention
- reduced wear ensures more efficiency

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Surgical cross-cut, particularly smooth-running, fast cutting of bone substance and tooth structure, **ZrN-coating as wear-protection coating**
total length: 45,0 mm

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Seperation of bone material, tooth root and teeth

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Fine 10-blade ball cutter, especially smooth-running

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Fine surgical cross-cut, particularly smooth-running

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Surgical cross-cut, particularly smooth-running, fast cutting of bone substance and tooth structure, total length: 65,0 mm

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
"Saw-toothing" with 4 blades for fast cutting
of bone substance, Total length: 70,0 mm

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
More instruments for surgery can be found in our MKG-catalogue.

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Pilot drill are used for the initial determination of the drilling depth. Afterwards the reaming of the drilling of the implant site is made by the final drill, for an optimal fit of the implant body.

- stainless special steel with a low cracking tendency
- sharp 2-blades-construction with big flutes for an easy chip removal
- s-shaped blade tip for easy penetrating
- relief-ground flutes for a low friction on the bone
- **laser-depth marking for accurate drilling (after each 2 mm), 2-blades after 8 mm, 3-blades after 6 mm drilling depth**

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Optimum speed: 800 min⁻¹ – 1.200 min⁻¹. Do not exceed the maximum of 2N.
Always use with spray, 50ml/min.

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For removal of the bone cover during the root tip resection.

- by the maximum drilling depth of 5.0 mm and the reduced total length, those instruments are also suitable for the posterior region (premolars/molars)
- their handiness and low drilling depth allow an optimal visibility on the work surface and reduce consequently the risk of injury
- these instruments are made of hardened stainless steel, for an easy identifying, the external and internal diameters are lasered on the shank
- easy grinding of the bone cylinders by the milled down windows on the instrument

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Recommendations for Trepan bur:

to use with contra-angle handpiece (minimum speed reduction 10:1) and permanent, sterile cooling.

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For the preparation of bone cylinders for bone harvesting.

- by the maximum drilling depth of 8 mm and the reduced total length, these instruments allow a better visibility on the work-surface and reduce consequently the risk of injury
- these instruments are made of hardened stainless steel, for an easy identifying, the external and internal diameters are lasered on the shank
- easy grinding of the bone cylinders by the milled down windows on the instrument

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For explanting.

- the maximum drilling depth allows a safe explanting of the commercially available implants up to the size of the diameter of 4.8 mm and a maximum total length of 18.0 mm
- the depth markings by laser allow an accurate working
- the milled down windows allow a good visibility on the implant body and support the cooling for an atraumatic explantation
- the instruments are made of hardened stainless steel, for an easy identifying, the external and internal diameters are lasered on the shank

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acurata® a-diamonds are manufactured according to state-of-the-art diamond bonding technology. The bonding layer is completely flat without any convexities. Each diamond grit is flatly embedded and cuts smoothly into the tooth substance. Electropolished cooling drainages allow the cooling fluid to circulate and provide easy and rapid removal of the grinding debris. The result is a constantly sharp instrument with long service life. Rounded head shapes secure round cavity line angles and crown abutments. Different grit sizes allow tailor-made roughness heights on the prepared surfaces.

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The cutting edges of the diamond grits are freely exposed. The lowered bonding layer holds a large amount of grinding debris. Each grit cuts smoothly into the tooth substance. **This is the acurata a-diamond.**

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A comparable instrument: Deeply embedded grits lead to low abrasive capacity. Dentine clogs the instrument. The resulting friction heats up the tooth and endangers the pulp.

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The Order Number describes the type of:

- material (working part)
- shank
- head shape
- grit size
- head diameter

Shape	250		Material	Shank	Shape	Grit	ISO Ø
Grit	 	=	806	314	250	534	018
ISO Ø	018		Diamond	FG	Flame	coarse	1,8 mm

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Ultra fine: 9 µm, Super finishing of composites / compomeres and ceramics.

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Extra fine: 25 µm, Finishing of composites, ceramics and filling materials.

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Fine: 40–54 µm, Finishing of prepared tooth surfaces and restorations.

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Medium: 90–126 µm, Standard applications cavity / crown / in- onlay preparation.

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Coarse: 154 µm, Rapid reduction of tooth substance and removal of old fillings.

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Super coarse: 180–200 µm, Coarse pregrinding / removal of old fillings.

Separation of none-precious alloys crowns and bridges.

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
- gracile working parts, slender necks of the instruments for an optimum sight on the working surface
- ideal in combination with magnifying glasses or microscopes
- the ideal grit size of 75 μm allows a precise preparation with smooth prepared surfaces at the same time
- minimal invasive shaping and working on the cavity
- maximum conservation of the healthy tooth substances
- practice-oriented selection of shapes

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acurata depth marker “type 102” is suitable for all types of crown- and bridge preparations. With the defined penetration depth in 4 lengths (0,5 mm – 2,0 mm) you can fix the exact preparation depth. Ideal for full ceramic restorations of crowns, bridges or veneers.

Tip: The length 2 mm is recommended for the occlusal full ceramic preparation.

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Particularities:

- Preparation of crowns, bridges
- defined depth, coating removal
- uniform wall thickness
- laser marked, no likelihood of confusion
- durable
- high cost efficient
- ideal for turbine, sufficient cooling (50 ml/min⁻¹)

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Prepares smooth and cuts extremely effective.

- high removal rate thanks to a special diamond coating
- coated and uncoated segments alternate with each other
- minimizes resonance vibration ("tuning fork effect")

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diamonds are gold coated. This feature allows them to dissipate
heat – even when large amounts of substance are cut – and it makes
the instruments absolutely corrosion-proof. **acurata gold** diamonds
are available in two grit sizes: 107–126 microns for smooth prepara-
tions and 50 microns for finishing. These instruments are ideal
for technically preparations in ceramic-fused-to-gold prosthetics.

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Facilitated lingual access in the lower jaw because of the extended shank. Figure: Removal of a lingual retainer in the lower jaw, coarse grit/green. 806 314 277LN 534 023

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Polishing of composite-restorations or removal of composite residues after a previous removal of a retainer with red grit. 806 314 277LN 514 023

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Extra-fine polishing of a composite-restoration with optimal sight due to the long slender shank. 806 314 277LN 504 023

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Collision with neighboring teeth and limited visibility during preparation or polishing with an egg with common (normal) shank length. 806 314 277 524 023

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Also facilitated preparation or treatment and polishing in the palatinal upper jaw front by the long slender shank. 806 314 277LN 514 023

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
Finishing of a palatinal crown preparation under optimal visibility (yellow grit). 806 314 277LN 504 023

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Good overview during a palatinal crown preparation thanks to a smooth running extended shank.
806 314 277LN 534 023

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Finishing of spots, which are difficult to reach thanks to the distally inclined bur head with excellent visibility. Wide range of applications, various uses, two different shapes in all available grits. 806 314 277LN 504 023

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a-diamond FG short, for the dentistry and patients with tight oral cavity and restricted mouth opening

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Recommendations for use: urbine, red contra-angle with mini-head low contact
pressure 0,2–2 N, cooling 50 ml/min.

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acurata has developed an ideal polish for each material used in dentistry. For amalgam, composites, compomeres, precious metals, ceramic-fused-to metal restorations with non-precious metals, silicates and feldspar ceramics, sintered oxide- and zircon oxide ceramics, resins, titanium. Polishing agents for prophylactic treatment lighten the teeth and smooth the enamel.

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Silicon carbide, green · ceramic bond, grit size medium for ceramic materials, metal alloys and composite

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Rough and rapid pre-polishing of all ceramics and composites. Rough and rapid removal of excess material. **Advantage for the user:** Very fast, efficient and optimal results by congruent form! Combinable with all acurata DIA-polishers, thus also including Pinky (composite) and the ceramic variants!

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For the polishing of Zircon oxide, Lithium-disilicate ceramics, silicate ceramics

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Smooth by a special bond, easily recognisable by the orange ring, optimum speed 5.000–10.000 min⁻¹.

We recommend previously finishing with the tungsten carbide finisher FQ (p. 22).

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Universally for all modern composites and compomer materials,

Artificial resin polishers interspersed with diamond grain – perfect result in two steps, don't smears – optimum speed 5.000–10.000 min⁻¹

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TURBO SHINE – DIAMOND POLISHING WHEEL WITH LAMELLAS DESIGN

- only three polishing steps to the perfect, prompt result on all ceramics
- diamond interspersed
- flexibly adapts to nearly every surface structure
- easy to use, work without pressure
- suitable for multiple uses and long lasting
- applicable with and without spray

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Prophy Color Brush Mini: For precise cleaning of surfaces and spots which are hard to reach e.g. multibracket appliances. Cleaning of fissures before sealing and lots of other applications.

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With integrated prophylaxis polishing paste. acurata Pasteless Prophy reaches a maximum cleaning and polishing power with minimum abrasion of the tooth substance

- latex free, with RA-shank
- high cleaning power for an optimum removal of stains
- reduced plaque formation, multiple use is possible

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New construction, napped structured surface, white plastic shank

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New shape, lamellar structured surface, flame shape

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ACURATA DISCOVERY ENDO-BUR – EXTRA LONG

- long slender neck
- extra long version, all-over length 34 mm
- to expand the root canal entrance
- for the Isthmus-preparation
- lasermarking 6 mm, 9 mm and 12 mm for the depth estimate
- color-coded for a secure recognition of the sizes, individual processing
- one-piece tungsten carbide construction – for an optimal stability
- max. speed: 20.000 min⁻¹, opt speed: 1.500 – 2.000 min⁻¹

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Longneck-Endo-Bur for the Isthmus-preparation.

- extra long and slender neck
- good sight on the working surface, safe in application
- sharp toothing
- clear identification by the green colour ring

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Recommendations: Optimum Speed: 1.000–2.000 min⁻¹, with light contact pressure. Do not lever or wedge the instrument.

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Recommendations for Use and Safe Operation of Root Canal Instruments:

- K-Reamer: push and rotary motion clockwise up to max. 90°
- K-File: filing motion clockwise up to max. 45°
- Hedstroem: push and pull (with or without rotary) motion up to max. 90°
- Usage in ascending order of size (do not skip any size)
- high flexible Nickel-Titanium instruments are especially for heavy curved root canals. The file has not to be bent before use
- protect the instruments against aspiration (Rubber dam – Safety chain – Safety cord)
- clean the instruments properly before sterilization
- Sterilization in autoclave may corrode steel instruments
- do not use bent or damaged instruments

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Color of stops indicates the length of working part

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
For safekeeping of rotary dental instruments in the workplace, during cleaning, disinfection and sterilisation. Safe fixing of the instruments by retainer clips. Universal elastic container holds FG, RA and HP.

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
Container for 8 FG or RA, holds instruments up to an all-over-length of 30 mm.

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
Container for 6 HP, holds instruments up to an all-over length of 40 mm.

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Container for 6 HP, holds instruments up to an all-over length of 50 mm.

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Container for 6 HP, holds instruments up to an all-over length of 60 mm.

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Recommendations:

- use only all known decontamination products for rotary instruments
- cleaning in ultrasonic bath is possible
- sterilisable up to 134 °C
- do not use hot air sterilizer and chemiclave because of the polymer parts
- put the bur stand into the cleaning device, so that the instruments are reached by the spray

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The "small range for laboratory" helps to finish immediately the necessary modifications.

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General Information for Use and Safety

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acurata instruments are intended for the dental medical field and must only be used by dentists or respective experts who, due their training and experience, are familiar with the use of the products. The instruments are supplied non-sterile and have to be prepared prior to first use (refer to Recommended Hygiene Procedures <http://www.acurata.de/en/service/recommended-hygiene-procedures>).

INFORMATION FOR USE FOR DENTAL ROTARY INSTRUMENTS

Use only technical and hygienic correct drives with a stable ball bearing! Carefully clamp the instruments as deep as possible and assure that the desired speed is reached prior to usage. Please note the maximum allowed engine speed as stated on the packaging label. Avoid levering or canting, otherwise there is risk of fracture. Pay attention that there is sufficient cooling (50 ml/min) otherwise there is risk of thermal damage. FG-instruments with a total length > 19 mm and instruments with a diameter of the work element > 2 mm (ISO 020) need extra cooling. Separate blunt, twisted or damaged instruments immediately and do not use them again. Improper usage increases risk, higher wear and inferior work results!

Contact pressures: The contact pressures should be light, do not exceed the maximum of 2N. Excessive contact pressures can cause damage to the instruments. In extreme cases the instrument can fracture! At the same time a higher heat development is generated. Overheating may harm the dental pulp, instruments wear faster and surfaces get coarser.

Remarks: Polishers and instruments with long working part and respectively long, slender necks (e.g. instruments for micro preparation) or for specific application areas (e.g. crown separator, long surgical cutter): exceeding the maximum allowed speeds can cause oscillations (vibrating) or displacement (whippings). This can result in fracture or heavy damage. Non observance causes an increased security risk! Therefore pay attention to the maximum speed printed on the packaging labels · The respective maximum speed is a theoretical value given on the basis of the geometry of the instruments. The choice of the working speed is determined by the material subsequently to be used, the specific indication, the type of handpiece used and the contact pressures and is at user's discretion. Please ensure a dry, dust-free and sun-protected storage at appropriate room temperature.

Maximum speeds – table of reference values

FG instruments (ISO 314) ISO Ø: 005–016 max. 450.000 min⁻¹ • 018–021 max. 300.000 min⁻¹ • 023–031 max. 160.000 min⁻¹ • 033–055 max. 120.000 min⁻¹

HP and RA instruments (ISO 104 / 204) ISO Ø: 003–023 max. 50.000 min⁻¹ • 025–040 max. 40.000 min⁻¹ • 045–060 max. 30.000 min⁻¹ • 060–070 max. 25.000 min⁻¹ • 070–250 max. 25.000 min⁻¹ • 300–380 max. 20.000 min⁻¹ • 450 max. 15.000 min⁻¹

Polishing instruments all ISO Ø: max. 5.000–10.000 min⁻¹

The information on the label on the packaging applies. Optimum speeds for all instruments depend on the material to be processed. As a general rule 40 – 50% of the maximum instruction is recommended!

RECOMMENDATIONS FOR USE FOR ROOT CANAL INSTRUMENTS

K-Reamer: push and rotary motion clockwise up to max. 90°. K-File: filing motion clockwise up to max. 45°. Hedstroem: push and pull (with or without rotary) motion up to max. 90°. Usage in ascending order of size (do not skip any size), High flexible Nickel-Titanium instruments are especially for heavy curved root canals. The file has not to be bent before use, Protect the instruments against aspiration (rubber dam – safety chain – safety cord), Silicone stoppers have to be removed prior to reprocessing root canal instruments, Clean the instruments properly before sterilization, Sterilization in autoclave may corrode steel instruments, Do not use bent or damaged instruments.

Cleaning, disinfection, sterilization and storage of medical products – refer to our Recommended Hygiene Procedures, also available on our website:
www.acurata.de/en-US/products/hygiene-recommendations.

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RECOMMENDED HYGIENE PROCEDURES FOR DENTISTRY

rev. 07/15

Manufacturer's information on the reprocessing of medical products semi critical B and critical B according to the RKI recommendation "Prevention of infections in dentistry hygiene requirements", the KRINKO guideline "Requirements on the hygiene for reprocessing of medical products" and the DIN EN ISO 17664.

Manufacturer: acurata GmbH & Co. KG · Schulstraße 25 · 94169 Thurmansbang · Tel.: +49(0)8504-91 17 0 · Fax.: +49(0)8504-91 17 90 · info@acurata.de · www.acurata.de

Products: The present manufacturer's information applies to all instruments supplied by acurata GmbH & Co. KG that are used for surgical, periodontal or endodontic treatments. These are rotating tungsten carbide, diamond and polishing instruments as well as oscillating instruments made of stainless steel or nickel-titanium. The product range includes exclusively instruments delivered in non-sterile condition. These have to be prepared prior to first use (begin with step 1) and prior to any further use (begin with step 2).

Limited number of reprocessing cycles: The end of a product's service life is determined by damage and wear caused by use. Frequent reprocessing does not affect the performance of these instruments. Disposable products (marked  on the label) shall not be reused or reprocessed.

Basic note: Observe the legal provisions concerning the reprocessing of medical products valid in your country (e.g. in Germany www.rki.de). The manufacturer assures that the reprocessing methods detailed below are appropriate for the reprocessing of the mentioned groups of instruments for their reutilization in regard of their intended use. The operator is responsible that the applied methods of reprocessing, based on his risk assessment with the used equipment, material, process parameters and staff achieve the required results for the intended use. To guarantee this, routine controls of the validated mechanical and/or manual reprocessing methods are necessary. Any deviation from the validated method below detailed must be checked and released by the operator to ensure effectiveness and to avoid possible adverse consequences. **The hygiene procedure recommendation is available on our website in its current version:**

<https://www.acurata.de/en-US/products/hygiene-recommendations>

1 Preparation incl. storage and transportation

For first use preparation begin with step 2. Place instruments immediately or at the latest one hour after use on a patient, in a cleaning/disinfection tank filled with a suitable detergent/disinfectant (non-fixing/ aldehyde-free, e.g. BIB forte eco). Set up the drill bath acc. to manufacturer's instructions; for BIB forte eco mix concentrate with water - first the water, then add concentrate. Cover the tank. Pay attention to the application time (e.g. BIB forte eco 0,5% 60 min.). The transport of the instruments to the place of preparation should be made in a contamination protected cleaning/disinfection tank.

2 Cleaning and disinfection

According to the directive of the Robert Koch Institute (RKI) and Commission for Hospital Hygiene and Infection Prevention (KRINKO) it is preferential that the preparation of semi-critical B products is carried out mechanically; critical B products shall in any case be prepared mechanically. For products with long, tight lumina or cavities the cleaning has to be performed mechanical. For root canal instruments silicone stoppers have to be removed prior to reprocessing.

Mechanical cleaning – validated method

Equipment: Cleaning brush, mechanical washer/disinfector acc. to EN ISO 15883 (RDG) (e.g. Miele with Vario TD program), detergent (e.g. 0,5% cleaner Neodisher mediclean), acurata bur stand made of stainless steel.

Method: Remove instruments from instrument stand / cleaning-disinfection-tank immediately before mechanical reprocessing and brush off all visible contamination in a cold water bath. Place the instruments into the opened bur stand. Start mechanical cleaning according to the instructions of the manufacturer of the machinery and of the detergent. The following process is validated: Program Vario TD: 2 min. precleaning, 5 min cleaning at 55°C with detergent, 3 min. neutralizing, 2 min. intermediate rinsing, final rinsing with appropriate (VE-) water 5 min. at >90°C.

Thermal disinfection in a validated washer/disinfector

Perform the mechanical cleaning (see above, e.g. Miele washer/disinfector with Vario TD program) incl. thermal disinfection with the instruments fixed in a bur holder. The manufacturer's instructions for the device must be observed. For validated washer/ disinfectors the disinfection is demonstrably assured. acurata products are thermostable up to 134°C.

Manual cleaning and disinfection – standardized method

Equipment: Cleaning brush (e.g. synthetic brush, sterilizable), ultrasonic bath, detergent & disinfectant with approved efficiency for dental instruments (e.g. BIB forte eco, Alpro Medical), bur stand for rotating oscillating instruments (e.g. acurata bur stand made of stainless steel); the manufacturers' instructions must be observed.

Method: Remove instruments from bur stand / cleaning-disinfection tank immediately before manual cleaning. Brush off all visible contamination in a cold water bath. Rinse the instrument and the bur stand under running water. Put the instrument into a suitable strainer element and place it into the ultrasonic unit filled with detergent & disinfectant. Perform cleaning and disinfection according to the instructions of the manufacturers of the ultrasonic bath and the detergent and disinfectant; e.g. BIB forte eco 3% - 10 min. at 55°C tested to EN 14476. After the application time rinse the instrument thoroughly with appropriate water (e.g. VE-water). Dry the instruments preferably with medical compressed-air. According to KRINKO the manual reprocessing is finalized by a thermal disinfection in a steam sterilizer. Follow manufacturer specifications.

Visual examination with a suitable enlarger to ensure that the instrument is clean and undamaged (an enlargement of 8x -10x is recommended). If after reprocessing still residues of contamination are visible, repeat the cleaning and disinfecting process until no visible contamination is left. Instruments showing defects are to be discarded immediately (e.g. missing diamond coating, blunt and chipped blades, deformations, corroded surfaces or non-removable residual contamination).

3 Final reprocessing steps – sterile packaging and sterilization – Medical products critical B – validated method with moist heat:

Equipment: Steam sterilizer Co. MMM Selectomat HP; acurata bur stand made of stainless steel, transparent sterilization bag (Steriking o. VP Stericlin), sealed seam device Co. Hawo

Packaging: Prior to sterilization place the instruments in the bur holder and pack them altogether doubly in sterilization bags and weld them with a sealed seam device. The instruments must be protected. For the packing an appropriate standardized method has to be applied.

Sterilization: An effective steam sterilization of the packaged instruments is proven successfully in the pre-vacuum steam sterilization method with the following minimal parameters: 3 pre-vacuum phases, 132°C sterilization temperature, holding time 3 min. (full cycle), drying time 10min. Follow the instructions of the device manufacturer. Note: The products are not suitable for a sterilization in a hot-air sterilizer or chemiclave.

4 Transport and storage During transport and storage, the reprocessed products must be protected from recontamination. Further the packed sterile goods must be protected also from dust and moisture.

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HIGH PERFORMANCE DENTAL INSTRUMENTS

to the highest technological level.

Dentists and dental laboratories now value our practically-oriented, needs-oriented products all over the world. We also listen and ask questions – at the same level as our customers. It is in this context that the highly focused, dedicated perfectionists at acurata manufacture absolutely reliable precision instruments which enable ultra-accurate work in the laboratory and in the dental practice. Reliability in the context of a sensitive collaboration between all of our talents – at the technical and human level.

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