

Zirkonoxid Schleifkörper

FG-Schaft, nur mit Spraykühlung verwenden.

acurata® hat mit dem „Zirkon-a-diamant“ FG-Instrumente entwickelt, deren feine Diamantkörnung in einer Spezialbindung exakt auf die Bearbeitung von gesinterter Zirkonoxidkeramik abgestimmt ist. Mikrorisse und thermische Überlastung werden vermieden. „Zirkon-a-diamant“ ist von ZrO₂-Anwendern erprobt.

Zircon oxide abrasives

FG shank, to be used only with spray cooling.

acurata designed the "zircon-a-diamond" FG instruments with a special fine diamond bonding especially tuned for sintered zircon oxide ceramics, to avoid micro cracks and overheating. The "zircon-a-diamond" is proved by ZrO₂-expert users.



Diamantkörnungen Diamond grit sizes
526 mittel medium 516 fein fine 506 extrafein extra fine 484 ultrafein ultra fine

Nur mit Spraykühlung verwenden!
Empfohlene Drehzahl: 160.000 min⁻¹ / Spray: 50 ml/min.

To be used only with spray cooling!
Recommended speed: 160.000 min⁻¹ / spray: 50 ml/min.

Allgemeine Bearbeitung

General application

Birne Pear



806 314 238 526 016

Flamme Flame



806 314 257 526 023



806 314 257 516 023

Granate Grenade



806 314 274 526 020



806 314 274 516 020

Konus Cone

Zylinder, rund Cylinder, round



806 314 142 526 017



806 314 142 516 017



806 314 142 506 017



806 314 142 484 017



806 314 142 526 023



806 314 142 516 023



806 314 142 506 023



806 314 142 484 023

23,5 mm

Konus Cone

Konus, rund Cone, round



806 314 199 526 020



806 314 199 516 020



806 314 199 506 020



806 314 199 484 020

23,5 mm

Konus Cone

Konus, rund Cone, round




806 314 199 526 023



806 314 199 516 023



806 314 199 506 023



806 314 199 484 023

23,5 mm



Zirkonoxid Schleifkörper

FG-Schaft, nur mit Spraykühlung verwenden.
 acurata® hat mit dem „Zirkon-a-diamant“ FG-Instrumente entwickelt, deren feine Diamantkörnung in einer Spezialbindung exakt auf die Bearbeitung von gesinterter Zirkonoxidkeramik abgestimmt ist. Mikrorisse und thermische Überlastung werden vermieden. „Zirkon-a-diamant“ ist von ZrO₂-Anwendern erprobt.

Zircon oxide abrasives

FG shank, to be used only with spray cooling.
 acurata designed the "zircon-a-diamond" FG instruments with a special fine diamond bonding especially tuned for sintered zircon oxide ceramics, to avoid micro cracks and overheating. The "zircon-a-diamond" is proved by ZrO₂-expert users.

Innenlumen Inner lumen

Zylinder, rund cylinder, round



806 314 L 137 526 016

Ei egg



806 314 L 277 526 014



806 314 L 277 516 014



806 314 L 277 506 014

26 mm

Kugel ball



806 314 697 526 012



806 314 697 526 016

24 mm

Allgemeine Bearbeitung General application

Zylinder, rund cylinder, round



806 314 140 526 012



806 314 142 526 014



806 314 142 516 014

Konus, rund cone, round



806 314 198 526 023



806 314 198 516 023



806 314 199 526 012



806 314 199 516 012



806 314 199 526 018



806 314 199 516 018

Diamantkörnungen Diamond grit sizes
 526 mittel medium 516 fein fine 506 extrafein extra fine 484 ultrafein ultra fine

Nur mit Spraykühlung verwenden!
 Empfohlene Drehzahl: 160.000 min⁻¹ / Spray: 50 ml/min.

To be used only with spray cooling!
 Recommended speed: 160.000 min⁻¹ / spray: 50 ml/min.